

# BMP™21-LAB Label Printer

*Tough on the outside, Genius on the inside*



Combining a tough exterior with smart printing capabilities, the BMP21-LAB hand-held label printer can take a hit and still perform. **Laboratory symbols** & vial/tube auto-sizing, plus drop-lock-and-go cartridges, plus a rechargeable Lithium-ion battery, plus a two year warranty adds up to a powerful and affordable tool for quick and efficient label making in any type of laboratory. **High quality label materials** designed specifically for the laboratory environment create durable, professional looking labels for your laboratory samples and beyond.

**Easy to use, easy to buy...you'll never use a pen again!**

## Designed for Laboratory

- Materials available for applications such as cryo vials and slides
- Greek symbols commonly used in laboratory, intuitive menu structure for symbols
- Vial label function automatically cuts the tape according to the size of the vial
- Automated serialisation
- Time & Date stamp

## Easy to use

- 3 line LCD display provides clear visibility
- ABC keypad, graphics library, hot keys and menu functions for fast label creation
- Long life rechargeable Lithium-ion battery\*
- Smart Cell Technology for automatic label setup
- Automatic label formatting for commonly used label types in laboratory
- Magnet\* for more convenient use and storage on smaller laboratory benches and a lanyard for hands-free operation
- Ergonomic, centre balanced printer with grab and go ridges grip for easy handling

## Exceptional material variety and durability

- Laboratory specific materials designed to withstand harsh laboratory environments.
- High quality, long lasting labels that stick and stay stuck
- Prints durable true-sized 6 – 19 mm wide identification to maximise amount of data
- Adjust the length of the labels using the high performance continuous material
- Fulfill your requirements with up to 7 different industry specified materials that last

**2** YEAR WARRANTY

\*Accessories sold separately. Battery included in the Starter Kit only.

## BMP™21-LAB Label Printer

Art. No.	Description
139538	BMP21-LAB Printer Blister Pack

**Blister package contains:** Printer + Lanyard + Quick Start Guide + M21-750-7425 label cartridge



**2 YEAR WARRANTY**

## BMP21-LAB Starter Kits

Art. No.	Description
134474	BMP21-LAB Printer Kit EU
134475	BMP21-LAB Printer Kit UK
134476	BMP21-LAB Printer Kit US

**Printer Kit contains:** Printer + Li-ion Battery + Power Cord + Lanyard + Quick Start Guide + M21-750-7425

## Accessories

Art. No.	Description
110416	AC Adaptor - Europe 220V
110417	AC Adaptor - UK 240V
110937	AC Adaptor - US 110V
139540	Li-ion battery Pack for BMP21-PLUS; BMP21-LAB
110891	Magnet Only Accessory
110892	Multifunctional Accessory

## Clear polyester labels for labelling of bottles, flasks and dishes



Order Ref.	Colour	Width (mm)	Length (m)
<b>B-430: Polyester</b>			
M21-375-430	Black on Clear	9.53	6.40
M21-500-430	Black on Clear	12.70	6.40
M21-750-430	Black on Clear	19.05	6.40

## Polypropylene labels for vials, tubes, slides and well plates

- Survives exposure to liquid nitrogen, freezer, water bath and autoclave.



Order Ref.	Colour	Width (mm)	Length (m)
<b>B-7425: Polypropylene</b>			
M21-375-7425	Black on White	9.53	6.40
M21-500-7425	Black on White	12.70	6.40
M21-750-7425	Black on White	19.05	6.40

## Ultra thin all white polyester labels for vials, tubes and conical tubes



Order Ref.	Colour	Width (mm)	Length (m)
<b>B-461: Self-laminating Polyester</b>			
M21-375-461	Black on White	9.53	6.40
M21-500-461	Black on White	12.70	6.40
M21-750-461	Black on White	19.05	6.40

## Vinyl labels for general purpose labelling

- Adheres to irregular, rough and curved surfaces
- Recommended for equipment, shelves, cabinets and general laboratory identification



Order Ref.	Colour	Width (mm)	Length (m)
<b>B-595: Vinyl</b>			
M21-375-595-GY	Black on Grey	9.53	6.40
M21-500-595-GY	Black on Grey	12.70	6.40
M21-750-595-GY	Black on Grey	19.05	6.40
M21-375-595-OR	Black on Orange	9.53	6.40
M21-500-595-OR	Black on Orange	12.70	6.40
M21-750-595-OR	Black on Orange	19.05	6.40
M21-250-595-WT	Black on White	6.35	6.40
M21-375-595-WT	Black on White	9.53	6.40
M21-500-595-WT	Black on White	12.70	6.40
M21-750-595-WT	Black on White	19.05	6.40
M21-250-595-YL	Black on Yellow	6.35	6.40
M21-375-595-YL	Black on Yellow	9.53	6.40
M21-500-595-YL	Black on Yellow	12.70	6.40
M21-750-595-YL	Black on Yellow	19.05	6.40
M21-375-595-BK	White on Black	9.53	6.40
M21-500-595-BK	White on Black	12.70	6.40
M21-750-595-BK	White on Black	19.05	6.40
M21-375-595-BL	White on Blue	9.53	6.40
M21-500-595-BL	White on Blue	12.70	6.40
M21-750-595-BL	White on Blue	19.05	6.40
M21-375-595-BR	White on Brown	9.53	6.40
M21-500-595-BR	White on Brown	12.70	6.40
M21-750-595-BR	White on Brown	19.05	6.40
M21-375-595-GN	White on Green	9.53	6.40
M21-500-595-GN	White on Green	12.70	6.40
M21-750-595-GN	White on Green	19.05	6.40
M21-375-595-PL	White on Purple	9.53	6.40
M21-500-595-PL	White on Purple	12.70	6.40
M21-750-595-PL	White on Purple	19.05	6.40
M21-375-595-RD	White on Red	9.53	6.40
M21-500-595-RD	White on Red	12.70	6.40
M21-750-595-RD	White on Red	19.05	6.40

## Chemical-resistant polyester labels for slides

- Withstands DMSO, ethanol and short term xylene exposure



Order Ref.	Colour	Width (mm)	Length (m)
<b>B-488: Polyester</b>			
M21-375-488	Black on White	9.53	6.40
M21-500-488	Black on White	12.70	6.40
M21-750-488	Black on White	19.05	6.40

## Nylon labels for bags, bottles and tubes

- Easy to handle nylon cloth label for non-rigid and curved surfaces



Order Ref.	Colour	Width (mm)	Length (m)
<b>B-499: Nylon Cloth</b>			
M21-375-499	Black on White	9.53	4.90
M21-500-499	Black on White	12.70	4.90
M21-750-499	Black on White	19.05	4.90



For more information on these materials visit [www.bradyeurope.com](http://www.bradyeurope.com).

**BRADY**

WHEN PERFORMANCE MATTERS MOST™

### Brady EMEA

Africa	+27 11 704 3295	Middle East	+971 4881 2524
Benelux	+32 (0) 52 45 78 11	Norway	+47 70 13 40 00
Central & Eastern Europe	+421 2 3300 4800	Romania	+40 21 202 3032
Denmark	+45 66 14 44 00	Russia	+7 495 225 93 62
France	+33 (0) 3 20 76 94 48	Spain & Portugal	+34 900 902 993
Germany	+49 (0) 6103 7598 660	Sweden	+46 (0) 8 590 057 30
Hungary	+36 23 500 275	Turkey	+90 212 264 02 20
Italy	+39 02 26 00 00 22	UK & Ireland	+44 (0) 1295 228 288

Your Brady Distributor

