



## Lab Coats

### Fisherbrand™ Basic Protection Polypropylene Disposable Lab Coats

These light and breathable latex-free, 50gsm, spunbond polypropylene lab coats feature a standard collar and five snaps for a secure closure.

#### White – Knit Cuffs, Three Pocket

Cat. No.	Size	Quantity
17-100-881A	S	Pack of 10
17-100-882A	M	Pack of 10
17-100-883A	L	Pack of 10
17-100-884A	XL	Pack of 10
17-100-885A	2XL	Pack of 10
17-100-886A	3XL	Pack of 10
17-100-887A	5XL	Pack of 10

#### Medical Blue – Knit Cuffs, Three Pocket

Cat. No.	Size	Quantity
17-100-889A	S	Pack of 10
17-100-890A	M	Pack of 10
17-100-891A	L	Pack of 10
17-100-892A	XL	Pack of 10
17-100-893A	2XL	Pack of 10
17-100-894A	3XL	Pack of 10
17-100-895A	5XL	Pack of 10

#### White – Elastic Cuffs, No Pockets

Cat. No.	Size	Quantity
17-100-801	S	Case of 25
17-100-802	M	Case of 25
17-100-803	L	Case of 25
17-100-804	XL	Case of 25
17-100-805	2XL	Case of 25
17-100-806	3XL	Case of 25
17-100-807	5XL	Case of 25



## Fisherbrand™ General Protection Disposable SMS Lab Coats

These 50gsm, SMS, single-use lab coats provide high durability with splash and dust protection. Features include knit cuffs, standard collar, three pockets and a five-snap closure to ensure a secure fit.

### White

Cat. No.	Size	Quantity
17-100-809A	S	Pack of 10
17-100-810A	M	Pack of 10
17-100-811A	L	Pack of 10
17-100-812A	XL	Pack of 10
17-100-813A	2XL	Pack of 10
17-100-814A	3XL	Pack of 10
17-100-815A	5XL	Pack of 10

### Teal

Cat. No.	Size	Quantity
17-100-817A	S	Pack of 10
17-100-818A	M	Pack of 10
17-100-819A	L	Pack of 10
17-100-820A	XL	Pack of 10
17-100-821A	2XL	Pack of 10
17-100-822A	3XL	Pack of 10
17-100-823A	5XL	Pack of 10

### Purple

Cat. No.	Size	Quantity
17-100-833A	S	Pack of 10
17-100-834A	M	Pack of 10
17-100-835A	L	Pack of 10
17-100-836A	XL	Pack of 10
17-100-837A	2XL	Pack of 10
17-100-838A	3XL	Pack of 10
17-100-839A	5XL	Pack of 10

### Pink

Cat. No.	Size	Quantity
17-100-849	S	Pack of 10
17-100-850	M	Pack of 10
17-100-851	L	Pack of 10
17-100-852	XL	Pack of 10
17-100-853	2XL	Pack of 10
17-100-854	3XL	Pack of 10
17-100-855	5XL	Pack of 10

### Cranberry

Cat. No.	Size	Quantity
17-100-857A	S	Pack of 10
17-100-858A	M	Pack of 10
17-100-859A	L	Pack of 10
17-100-860A	XL	Pack of 10
17-100-861A	2XL	Pack of 10
17-100-862A	3XL	Pack of 10
17-100-863A	5XL	Pack of 10

### Blueberry

Cat. No.	Size	Quantity
17-100-865A	S	Pack of 10
17-100-866A	M	Pack of 10
17-100-867A	L	Pack of 10
17-100-868A	XL	Pack of 10
17-100-869A	2XL	Pack of 10
17-100-870A	3XL	Pack of 10
17-100-871A	5XL	Pack of 10

### Medical Blue

Cat. No.	Size	Quantity
17-100-873A	S	Pack of 10
17-100-874A	M	Pack of 10
17-100-875A	L	Pack of 10
17-100-876A	XL	Pack of 10
17-100-877A	2XL	Pack of 10
17-100-878A	3XL	Pack of 10
17-100-879A	5XL	Pack of 10

**Assorted** (Includes: 1 each of White, Medical Blue, Pink, Cranberry; 2 each of Teal, Purple, Blueberry)


Cat. No.	Size	Quantity
17-100-841	S	Pack of 10
17-100-842	M	Pack of 10
17-100-843	L	Pack of 10
17-100-844	XL	Pack of 10
17-100-845	2XL	Pack of 10
17-100-846	3XL	Pack of 10
17-100-847	5XL	Pack of 10



## Fisherbrand™ Disposable Polypropylene Lab Coats

These white, 50gsm, disposable polypropylene lab coats provide breathable protection with a turndown collar, elastic cuff and velcro closure.

Cat. No.	Description	Quantity
17-100-940	S	Pack of 10
17-100-941	M	Pack of 10
17-100-942	L	Pack of 10
17-100-943	XL	Pack of 10
17-100-944	2XL	Pack of 10
17-100-945	3XL	Pack of 10
17-100-946	5XL	Pack of 10



## Fisherbrand™ Anti-Static Disposable SMS Lab Coats

Single use anti-static coated, 40gsm, SMS lab coats prevent static discharge, static jump, and static decay. Features include knit cuffs, standard collar, three pockets and a five-snap closure to ensure a secure fit.

### White

Cat. No.	Description	Quantity
17-100-926	S	Pack of 10
17-100-927	M	Pack of 10
17-100-928	L	Pack of 10
17-100-929	XL	Pack of 10
17-100-930	2XL	Pack of 10
17-100-931	3XL	Pack of 10
17-100-932	5XL	Pack of 10

### Medical Blue

Cat. No.	Description	Quantity
17-100-933	S	Pack of 10
17-100-934	M	Pack of 10
17-100-935	L	Pack of 10
17-100-936	XL	Pack of 10
17-100-937	2XL	Pack of 10
17-100-938	3XL	Pack of 10
17-100-939	5XL	Pack of 10

## Accessories



## Fisherbrand™ Protective Disposable Facemasks

Single-use latex free facemasks with a meltblown polypropylene filter that is thermally bonded between two layers of nonwoven polypropylene.

Cat. No.	Description	Quantity
17-100-896	White; Ear Loops; 2-ply	Case of 500
17-100-897	Blue; Ear Loops; 3-ply	Case of 500
17-100-898	White; Ear Loops; 3-ply	Case of 500
17-100-899	Blue; Ties; 3-ply	Case of 500



## Fisherbrand™ Disposable Poly Latex Shoe and Boot Covers

These disposable, 6mil, blue shoe and boot covers are constructed of polyethylene and polypropylene with elastic (latex-free) ankles.

Cat. No.	Description	Quantity
17-100-905	X-Large; Boot cover	Case of 200
17-100-909	Universal; Shoe cover	Case of 300
17-100-910	X-Large; Shoe cover	Case of 300



### Fisherbrand™ Polyethylene Disposable Shoe Covers

Get the protection of polyethylene in a disposable 3mil blue shoe cover.

Cat. No.	Description	Quantity
19-170-902	Universal	Case of 300
19-170-903	X-Large	Case of 300



### Fisherbrand™ Protective Disposable Shoe Covers

Get maximum protection with these disposable, blue 37gsm polypropylene shoe covers, available with or without anti-skid soles.

Cat. No.	Description	Quantity
17-100-903	Universal; Anti-Skid	Case of 300
17-100-904	X-Large; Anti-Skid	Case of 300
17-100-915	Universal; Without Anti-Skid	Case of 300
17-100-916	X-Large; Without Anti-Skid	Case of 300



### Fisherbrand™ Protective Disposable Anti-Static Shoe Covers

Sewn-in conductive carbon strips in these 42gsm disposable anti static polyethylene polypropylene shoe covers help eliminate static buildup while providing maximum protection.

Cat. No.	Description	Quantity
17-100-917	Universal	Case of 300



### Fisherbrand™ Disposable Non-Skid Shoe Covers

Low linting shoe covers made of 58gsm polyethylene/polypropylene provide optimal skid resistance on wet and smooth surfaces. Featuring blue color with white soles.

Cat. No.	Description	Quantity
17-100-911	Universal	Case of 200
17-100-912	X-Large	Case of 200



### Fisherbrand™ Protective Disposable High Tack Shoe Covers

These 106gsm polyethylene polypropylene coated slip resistant shoe covers provide superior protection.

Cat. No.	Description	Quantity
17-100-913	Universal	Case of 300
17-100-914	X-Large	Case of 300



### Fisherbrand™ Protective Polyethylene Sleeve Covers

Disposable, 6mil blue sleeve covers with latex-free elastic.

Cat. No.	Description	Quantity
19-170-904	Universal; 18" Length	Pack of 100



### Fisherbrand™ Basic Protection Disposable Pleated Bouffant Caps

The latex free elastic headband of these spunbonded 10gsm polypropylene caps is wrapped in soft polyester thread for great comfort and a snug fit Available in blue and white.

Cat. No.	Description	Quantity
17-100-920	Blue; 21"; Small	Pack of 100
17-100-922	Blue; 24"; Large	Pack of 100
17-100-919	White; 21"; Small	Pack of 100
17-100-921	White; 24"; Large	Pack of 100
17-100-923	White; 28"; X-Large	Pack of 100



### Fisherbrand™ Disposable Polypropylene Bouffant Caps

Get the protection of polypropylene in a comfortable, fitted 14gsm, latex free bouffant. Available in blue and white.

Cat. No.	Description	Quantity
19-170-900	White; 21"; Small	Pack of 100
19-170-901	White; 24"; Large	Pack of 100
17-100-900	Blue; 21"; Small	Pack of 100
17-100-901	Blue; 24"; Large	Pack of 100



### Fisherbrand™ Basic Protection Disposable Beard Covers

Cover beards and moustaches with white 12gsm spunbond polypropylene and latex-free elastic.

Cat. No.	Description	Quantity
17-100-925	Universal Size	Pack of 100

